ASTOCOVER

Passive warming insulation

The ASTOCOVER can be used in pre-operative, peri-operative and post-operative situations, in intensive care or on the ward, as prophylaxis against the patient cooling down and as a treatment for hypothermia. With ASTOCOVER it is prevented that warmth is radiated from the patient into the area and it comes thus to no unnecessary calorific loss of the patient. Hence the patient is protected against cooling down caused by operation theatre ventilation and draught.

An ideal solution to every application can be found in the above-mentioned scenarios by selecting applied parts which can be combined.

In combination with our warming systems such as ASTOPAD SYSTEM (Patient Warmer), ASTOFLO and ASTOTHERM PLUS (Blood and Infusion Warmers) a complete conception is realized to prevent the patient from the negative side effects of an operation such as hypothermia and decubitus.



ASTOCOVER PAS100 Upper body cover

to passive warming insulation, 100 X 100 cm

Order no.: PAS 100



Cover can be washed and disinfected

for ASTOCOVER PAS100, 100 X 100 cm

Order no.: COV 40100



ASTOCOVER PAS150 Body cover

to passive warming insulation, 50 X 150 cm

Order no.: PAS 150



Order no.: COV 40150



ASTOCOVER PAS155 Arm - Chest cover

to passive warming insulation, 50 X 150 cm

Order no.: PAS 155



Order no.: COV 40155



ASTOCOVER PAS180 Body cover

to passive warming insulation, 80 X 180 cm

Order no.: PAS 180



Order no.: COV 40180



ASTOCOVER PAS158 Arm-Upper body cover

to passive warming insulation, 100 X 150 cm

Order no.: PAS 158





Material:

The cover of ASTOCOVER consists of a PU coated PES fabric.

Characteristics:

- breathe-active
- virus and bacteria impermeable
- blood, secretion and urine impermeable

- flexible and no formation of wrinkles
- Ökotex 100, II Certificate

Cleaning and Disinfection:

For cleaning and disinfection usual, alcohol-based cleaning agents and disinfectants may be used. The washable and disinfectable covers for ASTOCOVER (COV 40100/40150/40155) additionally may be washed at 95°C.



LER ELECTRONIC

Stihler Flectronic GmbH

Julius-Hölder-Straße 36 • D-70597 Stuttgart • Germany www.stihlerelectronic.de • E-Mail: info@stihlerelectronic.de national: Telefon 07 11- 72 06 70 • Telefax 07 11- 720 67 57

international: Telefon ++ 49- 711- 72 06 70 • Telefax ++ 49- 711- 720 67 57